Open-Source Syringe Pump

(developed by Andrey Samokhin in 2018)

Assembly Instructions

(Version 1.0a)

STL-files and "firmware" are available at http://www.mass-spec.ru/projects/diy/syringe_pump/eng/
Project files are distributed under CC BY-SA 3.0 license






Fig. 1. Electrical circuit.





Fig. 2. Electronic components.

Step 1.1 (if necessary, change the basic settings in the Arduino sketch)

<u>Step 1.2</u> (upload sketch to Arduino) Detailed instructions are available at this link.

Step 1.3 (solder wires and crimp connectors, refer to electrical circuit)

PART#	PART NAME	QTY
22	Arduino UNO	1
26	Active Buzzer (D=12mm)	1
37	Dupont Connector (Female)	6
38	Dupont Connector (Male)	3
39	Dupont Housing (3P)	2
40	Connector Wires	_



Step 1.4 (solder wires and crimp connectors, refer to electrical circuit)

PART#	PART NAME	QTY
11	Endstop Switch	2
38	Dupont Connector (Male)	3
39	Dupont Housing (3P)	1
40	Connector Wires	_




Step 1.5 (plug A4988 driver into the shield and set all DIP switches to ON)

PART#	PART NAME	QTY
24	Stepper Motor Driver (A4988)	1
25	Control Shield for A4988	1




<u>Step 1.6</u> (set output current limit on A4988) Detailed instructions are available at this link.

<u>Step 1.7</u> (connect all electronic components together; connect stepper motor and endstop switches; plug in 12 V power supply; make sure that pump works)




Step 1.8 (place active buzzer, Dupont housing, and shield with A4988 into the case)

PART#	PART NAME	QTY
19	Case (Base)	1
35	Self-tapping Screw 2x8	2




Step 1.9 (place Arduino UNO into the case)

PART #	PART NAME	QTY
35	Self-tapping Screw 2x8	2



Step 1.10 (plug LCD Keypad Shield into Arduino UNO)

- 1	<u> </u>		
	PART#	PART NAME	QTY
	23	LCD Keypad Shield	1



Step 1.11 (install buttons and top cover)

PART#	PART NAME	QTY
20	Case (Cover)	1
21	Buttons	1
36	Self-tapping Screw 2x10	4



